

K
a
v
i
n

S
c
i
e
n
t
i
f
i
c


P
r
o
d
u
c
t
s



Kavin Scientific Products

PRICE LIST

2015-2016





Kavin Scientific Products

PRICE LIST

CHENNAI

7A. Natarajan St., (Saidai Sambadhan St.)
Saidapet, Chennai-15.
Phone : 24714329, Tele Fax : 044-43302975
Mobile : 98400 54224, 99400 94224
Email : kavinscientific@yahoo.co.uk

MADURAI

322. South Marat Street.
South Gate, Madurai-625001
Mobile : 98404 54224
Email : kavinscientific.mdu@gmail.com



LABORATORY GLASSWARES

(MADE FROM BOROSILICATE GLASS)



This type of borosilicate glass contains approx. 81% silicon dioxide and approx. 13% boron trioxide is featured in high chemical durability, low coefficient of linear thermal expansion and therefore even to high resistance against temperature changes

Borosilicate glass represent an internationally defined glass type, suitable for manufacturing the majority of laboratory glass, devices and industrial apparatus, The borosilicate glass features : Excellent corrosion resistance, smooth unporous surface, perfect transparency, catalytically indifferent, odourless and tasteless,

CHEMICAL DURABILITY

Borosilicate glass resists practically all aggressive agent except Hydrofluoric acid, Phosphoric acid and Concentrated Hot alkalis which considerably affect the surfaces of the glass. Hydrofluoric acid damages glass surface mostly at low concentrations and it is thus applied to etching scales on volumetric glass. Phosphoric acid and alkalis damage glass only slightly at low temperatures and concentrations. Continual changing aciduous and alkalinous medium increases corrossions.


REVIEW OF THE MOST IMPORTANT PHYSICAL CHARACTERISTICS

Medium linear coefficient of thermal expansivity α	20-300° C	3,3 10 ⁻⁶ K ⁻¹
Water resistance	98° C	class 1
Softening point (Littleton's)	10 ^{7.6} PAS	825°C
Density	20°C	2.23 g cm ⁻³




**1001-Beakers:Low form,with Graduation & spout,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
5	20	25	10	35.00
10	25	32	10	37.00
25	34	45	10	39.00
50	40	55	10	40.00
100	50	73	10	42.00
150	55	80	10	48.00
250	65	98	6	55.00
400	80	110	6	89.00
500	85	120	6	95.00
600	90	125	6	125.00
1000	105	145	6	190.00
2000	130	190	4	450.00



**1003-Beakers:Tall form with graduation & spout,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
50	38	70	12	42.00
100	48	80	12	45.00
150	54	95	12	46.00
250	60	120	6	70.00
400	70	130	6	92.00
500	75	140	6	95.00



2002 - Bottles, B.O.D with I/C Stopper, Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
125	55	120	12	235.00
300	70	150	12	325.00



**2004 - Bottles,Reagent Narrow Mouth with I/C Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	140.00
125	55	95	12	170.00
250	65	150	6	185.00
500	80	170	6	270.00
1000	100	200	2	420.00
2000	130	265	1	675.00



**2006 - Bottles,Reagent Narrow Mouth with I/C Stopper, AMBER
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	195.00
125	55	95	12	250.00
250	65	150	6	300.00
500	80	170	6	450.00
1000	100	200	2	670.00
2000	130	265	1	850.00





**2008-Bottles,Reagents Wide Mouth with I/C Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	170.00
125	55	95	12	185.00
250	65	150	6	235.00
500	80	170	6	410.00
1000	100	200	2	420.00
2000	130	265	1	700.00

**2010-Bottles,Reagents Wide Mouth with I/C Stopper, AMBER
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	45	95	12	225.00
125	55	95	12	247.00
250	65	150	6	320.00
500	80	170	6	575.00
1000	100	200	2	775.00
2000	130	265	1	820.00

**2012 - Bottles, Laboratory Graduated,With Screwcap and Pouring Ring,
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
60	36	74	10	120.00
50	46	91	10	135.00
100	56	105	10	225.00
250	70	143	10	250.00
500	86	181	10	290.00
1000	101	230	10	425.00
2000	136	267	6	735.00

**2014 - Bottles, Laboratory Graduated,With Screwcap and Pouring Ring, AMBER
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
25	36	74	10	175.00
50	46	91	10	196.00
100	56	105	10	260.00
250	70	143	10	425.00
500	86	181	10	520.00
1000	101	230	10	650.00
2000	136	267	6	1000.00

**2016 - Bottles,Specific Gravity with I/C Capillary bore Glass Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN**




Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
10	24	45	10	90.00
25	38	55	10	110.00
50	49	65	6	105.00



**2018 - Bottles, Specific Gravity with I/C Capillary bore Teflon Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
10	24	45	10	220.00
25	38	55	10	260.00
50	49	65	6	295.00



2020 - Washbottles Complete : Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
250	4	230.00
500	4	310.00
1000	2	395.00



**2022 - Bottles, Weighing Tall form with I/C Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
5	20	40	10	72.00
15	25	60	10	83.00
25	30	65	10	110.00
60	40	90	6	150.00



**2024 - Bottles, Weighing Squat form with I/C Stopper,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
20	50	35	10	230.00
40	65	40	10	274.00



**3007 - Burettes, Straight bore Glass Solid Stopcock : Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Grad.Civ	Tol+mm	Qty/pack	Price/Pc
10	0.05	0.05	10	170.00
25	0.10	0.10	10	180.00
50	0.10	0.10	10	185.00
100	0.10	0.10	4	200.00



**3009 - Burettes, Rotoflow Telfon Screwcock : Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Grad.Civ	Tol+mm	Qty/pack	Price/Pc
10	0.05	0.05	10	170.00
25	0.10	0.10	10	180.00
50	0.10	0.10	10	185.00
100	0.10	0.10	4	200.00



**3011 - Burettes, Straight bore Telfon Stopcock : Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Grad.Civ	Tol+mm	Qty/pack	Price/Pc
10	0.05	0.05	10	365.00
25	0.10	0.10	10	380.00
50	0.10	0.10	10	397.00
100	0.10	0.10	4	425.00





**3001 - Bottles, Straight bore Glass Stopcock : Class B - SCHELLBACH,
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Grad.Civ	Tol+mm	Qty/pack	Price/Pc
10	0.05	0.05	10	180.00
25	0.10	0.10	10	195.00
50	0.10	0.10	10	207.00
100	0.10	0.10	4	240.00

**3003 - Bottles, rotoflow Telfon Screw Cock : Class B - SCHELLBACH,
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Grad.Civ	Tol+mm	Qty/pack	Price/Pc
10	0.05	0.05	10	180.00
25	0.10	0.10	10	195.00
50	0.10	0.10	10	207.00
100	0.10	0.10	4	240.00

**3005 - Bottles, Straight bore Telfon Stopcock : Class B - SCHELLBACH,
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Grad.Civ	Tol+mm	Qty/pack	Price/Pc
10	0.05	0.05	10	475.00
25	0.10	0.10	10	495.00
50	0.10	0.10	10	545.00
100	0.10	0.10	4	690.00

**3013 - Burettes, Automatic, With Reservoir, Straight bore Glasskey Stopcock Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Qty/pack	Price/Pc
10	1	2300.00
25	1	2400.00
50	1	2650.00
100	1	2900.00

**4001 - Condensers, Air with I/C inner and outer Joint,
Made from Heat Resistant Borosilicate Glass, KAVIN**



App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
100	10/19	4	130.00
200	14/23	4	200.00
200	19/26	4	200.00
300	19/26	4	250.00
400	19/26	4	295.00

**4002 - Condensers, Liebig with I/C inner and outer Joint,
Made from Heat Resistant Borosilicate Glass, KAVIN**



App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
200	B19	1	230.00
300	B24	1	260.00
400	B24	1	310.00



**4005 - Condensers, Allihn with I/C inner and outer Joint,
Made from Heat Resistant Borosilicate Glass, KAVIN**

App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
200	B19	1	260.00
300	B24	1	290.00
400	B24	1	430.00



**4007 - Condensers, Graham Coil with I/C inner and outer Joint,
Made from Heat Resistant Borosilicate Glass, KAVIN**

App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
300	B24	1	460.00
400	B24	1	570.00
500	B24	1	650.00




**4009 - Condensers, Double Surface Davies with I/C inner and outer Joint,
Made from Heat Resistant Borosilicate Glass, KAVIN**

App.Length of Jacket mm	Joint Size	Qty/pack	Price/Pc
150	B19	1	490.00
200	B19	1	520.00
300	B19	1	550.00
300	B24	1	660.00
300	B29	1	750.00



5001 - Cylinder Nessler, Colour Comparison Made from Heat Resistant Borosilicate Glass

Cap.ml	Qty/pack	Price/Pc
50ml	10	68.00
100ml	10	95.00



**5003 - Cylinder Measuring Class B, Round base
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Grad.Civ	Tol $\pm$ mm	Qty/pack	Price/Pc
5	0.1	0.1	10	55.00
10	0.2	0.2	10	60.00
25	0.5	0.5	10	85.00
50	1.0	1.0	10	100.00
100	1.0	1.0	10	130.00
250	2.0	2.0	10	225.00
500	5.0	5.0	10	270.00
1000	10.0	10.0	4	371.00
2000	20.0	10.0	2	677.00





5005 - Cylinder Measuring Class B, Hexagonal base
Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Grad.Civ	Tol $\pm$ mm	Qty/pack	Price/Pc
5	0.1	0.1	10	130.00
10	0.2	0.2	10	140.00
25	0.5	0.5	10	160.00
50	1.0	1.0	10	180.00
100	1.0	1.0	10	215.00
250	2.0	2.0	10	380.00
500	5.0	5.0	10	510.00
1000	10.0	10.0	4	825.00
2000	20.0	10.0	2	1400.00

5007 - Cylinder Measuring Class B, Hexagonal base, Schellbach
Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Grad.Civ	Tol $\pm$ mm	Qty/pack	Price/Pc
5	0.1	0.1	10	145.00
10	0.2	0.2	10	160.00
25	0.5	0.5	10	183.00
50	1.0	1.0	10	225.00
100	1.0	1.0	10	245.00
250	2.0	2.0	10	440.00
500	5.0	5.0	10	575.00
1000	10.0	10.0	4	895.00
2000	20.0	10.0	2	1500.00



6001 - Dessiccators, Plain with Cover Knob Top KAVIN, Toughened Natural Glass

ID mm	Qty/pack	Price/Pc
150	1	2300.00
210	1	2550.00
240	1	3300.00
300	1	4100.00



6003 - Dessiccators, Vaccum with Cover Knob Top KAVIN, Toughened Natural Glass

ID mm	Qty/pack	Price/Pc
150	1	4200.00
210	1	4400.00
240	1	5800.00
300	1	8500.00



7001 - Dishes Petri, Made from Heat Resistant Borosilicate Glass, KAVIN

Approx O.D.xht mm	Qty/pack	Price/Pc
50x20	20	60.00
80x20	10	75.00
100x20	10	83.00
150x20	10	190.00



7003 - Dean & Stark Apparatus, Made from Heat Resistant Borosilicate Glass, KAVIN

Specefication	Qty/pack	Price/Pc
Receiver 10ml, R.B. Flask 1000 ml. Liebig Condenser,	1	1050.00



**7006 - Extraction Apparatus, Soxhlet Complete with allihn Condenser and F.B. Flasks
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Flask Cap.ml	Qty/pack	Price/Pc
60	150	1	800.00
100	250	1	1200.00
200	500	1	1625.00



**7007 - Extraction Apparatus, Spare Alliah Condenser
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Joint Bottom	Qty/pack	Price/Pc
For 60	34/35	1	250.00
For 100	40/38	1	350.00
For 200	50/42	1	650.00



**7008 - Extraction Apparatus, Spare Extractor
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	I/C Joint Top	I/C Joint Bottom	Qty/pack	Price/Pc
60	34/35	24/29	1	370.00
100	40/38	24/29	1	550.00
200	50/42	24/29	1	700.00



8001 - Flasks Conical, Narrow Mouth Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
10	34	60	20	45.00
25	42	70	20	46.00
50	51	82	20	47.00
100	64	105	10	55.00
150	74	120	10	65.00
250	84	140	6	78.00
500	104	175	6	120.00
1000	130	220	6	180.00
2000	163	285	1	430.00



**8002 - Flasks Conical, wide Mouth with Graduation
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
100	64	105	10	63.00
250	72	124	6	87.00
500	104	175	6	135.00
1000	131	210	6	210.00



8003 - Flasks Conical, with I/C Stopper Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	I/C Stopper	Qty/pack	Price/Pc
100	B19	10	125.00
100	B24	10	130.00
250	B24	6	140.00
250	B29	6	145.00
500	B24	4	225.00
500	B29	4	270.00
1000	B24	2	350.00
1000	B29	2	385.00
2000	B24	1	650.00
2000	B29	1	695.00
2000	B34	1	715.00





8004 - Flasks Iodine, with I/C Stopper Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	I/C Stopper	Qty/pack	Price/Pc
100	B24	10	110.00
250	B24	6	160.00
500	B24	4	230.00
1000	B24	2	380.00

**8005 - Flasks Flat Bottom, Narrow Neck Plain
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
25	42	90	10	48.00
50	51	95	10	50.00
100	64	110	10	60.00
150	75	137	10	74.00
250	85	145	10	75.00
500	105	178	6	115.00
1000	133	212	6	170.00
2000	167	230	2	333.00

**8006 - Flasks Round Bottom, Narrow Neck Plain
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
25	42	90	10	48.00
50	51	95	10	50.00
100	64	110	10	60.00
150	75	137	10	74.00
250	85	145	10	75.00
500	105	178	6	115.00
1000	133	212	6	170.00
2000	167	230	2	333.00

**9001 - Flasks, Flat Bottom, Centre Neck with I/C Joint
Made from Heat Resistant Borosilicate Glass, KAVIN**



Cap.ml	Joint Size	Qty/pack	Price/Pc
50	19/26	10	100.00
50	24/29	10	100.00
50	29/32	10	100.00
50	34/35	10	100.00
100	19/26	10	120.00
100	24/29	10	120.00
100	29/32	10	120.00
100	34/35	10	120.00
150	19/26	8	135.00
150	24/29	8	135.00
150	29/32	8	135.00
150	34/35	8	135.00
250	19/26	8	175.00
250	24/29	8	175.00
250	29/32	8	175.00
250	34/35	8	175.00
500	19/26	4	225.00
500	24/29	4	225.00
500	29/32	4	225.00
500	34/35	4	225.00



Cap.ml	Joint Size	Qty/pack	Price/Pc
1000	19/26	4	290.00
1000	24/29	4	290.00
1000	29/32	4	290.00
1000	34/35	4	290.00
2000	19/26	2	550.00
2000	24/29	2	550.00
2000	29/32	2	550.00
2000	34/35	2	550.00



**9002 - Flasks, Round Bottom, Centre Neck with I/C Joint
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Joint Size	Qty/pack	Price/Pc
50	19/26	10	100.00
50	24/29	10	100.00
50	29/32	10	100.00
50	34/35	10	100.00
100	19/26	10	120.00
100	24/29	10	120.00
100	29/32	10	120.00
100	34/35	10	120.00
150	19/26	8	135.00
150	24/29	8	135.00
150	29/32	8	135.00
150	34/35	8	135.00
250	19/26	8	175.00
250	24/29	8	175.00
250	29/32	8	175.00
250	34/35	8	175.00
500	19/26	4	225.00
500	24/29	4	225.00
500	29/32	4	225.00
500	34/35	4	225.00
1000	19/26	4	290.00
1000	24/29	4	290.00
1000	29/32	4	290.00
1000	34/35	4	290.00
2000	19/26	2	550.00
2000	24/29	2	550.00
2000	29/32	2	550.00
2000	34/35	2	550.00



**9003 - Flasks, Round Bottom, Two Necks, Centre Neck and One Side Angled Neck with
I/C Joint. Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Centre Neck	Centre Neck	Qty/pack	Price/Pc
100	24/29	14/23	5	175.00
100	24/29	19/26	5	175.00
250	29/32	14/23	4	220.00
250	34/35	19/26	4	220.00
500	19/26	14/23	2	250.00
500	24/29	19/26	2	250.00
1000	29/32	14/23	1	325.00
1000	34/35	19/26	1	325.00
2000	29/32	14/23	1	650.00
2000	34/35	19/26	1	650.00





9004 - Flasks, Round Bottom, Two Necks, Centre Neck and Parallel Neck with I/C Joint.
Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Centre Neck	Side Neck	Qty/pack	Price/Pc
100	24/29	14/23	2	175.00
100	24/29	19/26	2	175.00
250	29/32	14/23	2	220.00
250	34/35	19/26	2	220.00
500	19/26	14/23	2	250.00
500	24/29	19/26	2	250.00
1000	29/32	14/23	1	325.00
1000	34/35	19/26	1	325.00
2000	29/32	14/23	1	650.00
2000	34/35	19/26	1	650.00

9005 - Flasks, Round Bottom, Three Necks, Centre Neck and Two Angled Neck with I/C Joint
Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Centre Neck	Side Neck	Qty/pack	Price/Pc
100	14/23	14/23	2	229.00
100	19/26	14/23	2	240.00
100	24/29	14/23	2	240.00
100	24/29	19/26	2	245.00
250	19/26	14/23	2	290.00
250	19/26	19/26	2	290.00
250	24/29	14/23	2	290.00
250	24/29	19/26	2	290.00
250	29/32	19/26	2	375.00
500	19/26	14/23	2	375.00
500	19/26	19/26	2	375.00
500	24/29	14/23	2	375.00
500	24/29	19/26	2	375.00
500	29/32	14/23	2	390.00
500	29/32	19/26	2	390.00
1000	24/29	14/23	1	610.00
1000	24/29	19/26	1	610.00
1000	34/35	19/26	1	640.00
2000	24/29	14/23	1	764.00
2000	24/29	19/26	1	764.00
2000	34/35	19/26	1	764.00

9006 - Flasks, Filter Conical Shape with Glass Hose Connection
Made from Heat Resistant Borosilicate Glass, KAVIN



Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
100	64	105	10	131.00
250	85	150	6	220.00
500	105	184	4	290.00
1000	136	245	2	520.00
2000	168	295	1	1100.00



9007 - Flasks, Distilling with side arm Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
100	64	210	6	140.00
250	68	214	4	200.00
500	88	226	2	225.00
1000	110	250	1	340.00



9008 - Flasks, Kjeldahi Plain Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Approx O.D.	Ht mm	Qty/pack	Price/Pc
300	97	335	6	200.00
500			6	245.00
800	110	375	6	325.00




9009 - Flasks, Kjeldahi with i/C Joint, Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Joint Size	Qty/pack	Price/Pc
300	B24	6	300.00
500	B24	6	335.00
800	B24	6	435.00



**1010 - Flasks, Volumetric with I/C Stopper Class B
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Tel±mm	Joint	Qty/pack	Price/Pc
5	0.05	B10	10	75.00
10	0.05	B10	10	80.00
20	0.08	B10	10	85.00
25	0.08	B10	10	92.00
50	0.12	B10	10	99.00
100	0.20	B14	10	125.00
200	0.30	B14	10	140.00
250	0.30	B14	6	145.00
500	0.50	B19	6	199.00
1000	0.80	B19	6	315.00



**1020 - Flasks, Volumetric with I/C Stopper Class B AMBER
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Tel±mm	Joint	Qty/pack	Price/Pc
5	0.05	B10	10	103.00
10	0.05	B10	10	123.00
20	0.08	B10	10	129.00
25	0.08	B10	10	133.00
50	0.12	B10	10	145.00
100	0.20	B14	10	215.00
200	0.30	B14	10	225.00
250	0.30	B14	6	245.00
500	0.50	B19	6	378.00
1000	0.80	B19	6	593.00



2100 - Funnels, 60 Long Stem, Made from Heat Resistant Borosilicate Glass, KAVIN

Dia.mm	Qty/pack	Price/Pc
50	6	42.00
65	6	53.00
75	6	67.00
100	2	97.00



**2110 - Funnels, 60 Short Stem, for Powder,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Dia.mm	Qty/pack	Price/Pc
125	4	159.00
150	4	199.00
200	2	388.00





2120 - Funnels, Separating Pear Shape with Rototeflon Teflon Screwcock & I/C Stopper
Made from Heat Resistant Borosilicate Glass, KAVIN

Dia.mm	Qty/pack	Price/Pc
60	4	160.00
125	4	205.00
250	4	249.00
500	2	380.00
1000	2	792.00



2130 - Funnels, Separating Pear Shape with Glass Stopcock & I/C Stopper
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
60	4	130.00
125	4	175.00
250	4	225.00
500	2	310.00
1000	2	475.00



2140 - Funnels, Separating Pear Shape with Teflon Stopcock & I/C Stopper
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
60	4	495.00
125	4	523.00
250	4	555.00
500	2	745.00
1000	2	1100.00



2150 - Funnels, Separating Globe Shape with Rototeflon Teflon Screwcock & I/C Stopper
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
60	4	160.00
125	4	205.00
250	4	249.00
500	2	380.00
1000	2	792.00



2160 - Funnels, Separating Globe Shape with Glass Stopcock & I/C Stopper
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
60	4	145.00
125	4	195.00
250	4	235.00
500	2	365.00
1000	2	725.00



2170 - Funnels, Separating Globe Shape with Teflon Stopcock & I/C Stopper
Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
60	4	495.00
125	4	523.00
250	4	525.00
500	2	745.00
1000	2	1100.00



**3110 - Filter Crucibles with Sintered Disc,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Porosity Grade	Qty/pack	Price/Pc
15	G1	10	95.00
15	G2	10	95.00
15	G3	10	95.00
15	G4	10	95.00
30	G1	10	130.00
30	G2	10	130.00
30	G3	10	130.00
30	G4	10	130.00
50	G1	10	175.00
50	G2	10	175.00
50	G3	10	175.00
50	G4	10	175.00



**3120 - Funnels Buchner with Sintered Disc,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Porosity Grade	Qty/pack	Price/Pc
35	G1	2	165.00
35	G2	2	165.00
35	G3	2	165.00
35	G4	2	165.00
80	G1	2	330.00
80	G2	2	330.00
80	G3	2	330.00
80	G4	2	330.00
200	G1	2	525.00
200	G2	2	525.00
200	G3	2	525.00
200	G4	2	525.00



**4110 - Pipettes, Graduated, (Serological Type) Graduated to Tip Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Grad. Interval mm	Tol+ml	Qty/pack	Price/Pc
1	0.01	0.01	10	45.00
1	0.10	0.01	10	45.00
2	0.02	0.02	10	48.00
2	0.10	0.02	10	48.00
5	0.05	0.05	10	55.00
5	0.10	0.10	10	55.00
10	0.10	0.10	10	65.00
25	0.10	0.20	10	75.00



**4120- Pipettes, Graduated, (Mohrs Type) Graduated to Tip Class B,
Made from Heat Resistant Borosilicate Glass, KAVIN**

Cap.ml	Grad. Interval mm	Tol+ml	Qty/pack	Price/Pc
1	0.01	0.01	10	45.00
1	0.10	0.01	10	45.00
2	0.02	0.02	10	48.00
2	0.10	0.02	10	48.00
5	0.05	0.05	10	55.00
5	0.10	0.10	10	55.00
10	0.10	0.10	10	65.00
25	0.10	0.20	10	75.00





4130 - Pipettes, Volumetric, Class B, Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Tol+ml	Qty/pack	Price/Pc
1	0.015	10	45.00
2	0.020	10	45.00
5	0.030	10	48.00
10	0.040	10	48.00
20	0.060	10	55.00
25	0.060	10	55.00
50	0.100	10	65.00
100	0.160	10	76.00



5110 - Culture Tubes, Round Bottom with PP Screw Cap, Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	17.00
15	100	21.00
30	100	27.00
60	100	48.00



5120 - Culture Tubes, Flat Bottom with PP Screw Cap, Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	17.00
15	100	21.00
30	100	28.00



5130 - Culture Tubes, Round Bottom with PP Screw Cap, AMBER Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	42.00
15	100	72.00
30	100	82.00
60	100	92.00



5140 - Culture Tubes, Flat Bottom with PP Screw Cap, AMBER Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
5	100	42.00
15	100	59.00
30	100	65.00



5150 - Centrifuge Tube : Conical Plain Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
15	100	20.00
50	100	60.00




5160 - Centrifuge Tube : Conical Graduated Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Qty/pack	Price/Pc
15	100	28.00
50	100	90.00



5170 - Blood Sugar Tubes Made from Heat Resistant Borosilicate Glass, KAVIN

Cap.ml	Graduation	Price/Pc
25	12.5 & 25	44.00
25	4.12.5 & 25	52.00



6110 - Solid Stoppers : Made from Heat Resistant Borosilicate Glass, KAVIN

Size	Qty/pack	Price/Pc
B12	25	40.00
B14	25	40.00
B19	25	45.00
B24	25	75.00
B29	25	95.00
B34	25	135.00
B40	25	190.00
B50	25	200.00
B55	25	210.00



7110 - Test Tubes with Rim Made from Heat Resistant Borosilicate Glass, KAVIN

Size	Qty/pack	Price/Pc
10X75	100	6.00
12X75	100	7.00
12X100	100	8.00
15X125	100	9.00
15X150	100	10.00
18X150	100	11.00
25X150	100	22.00
32X200	25	45.00



7120 - Test Tubes without Rim Made from Heat Resistant Borosilicate Glass, KAVIN

10X75mm	100	6.00
12X75mm	100	7.00
12X100mm	100	8.00
15X125mm	100	9.00
15X150mm	100	10.00
18X150mm	100	11.00
25X150mm	100	22.00
32X200mm	25	40.00



**8110 - Membrane filter Holder all glass Assembly
Made from Heat Resistant Borosilicate Glass, KAVIN**

Funnel for 47mm dia filter base with sintered disc & Buchner Flask Cap. 1000ml. Funnel 300ml with suitable clamp.	3500.00
--	---------





LIQUID HANDLING PRODUCTS - BIOPETTE +

MICROPIPETTES, Fixed Volume Single channel, Fully Autoclavable, BIOPETTE Made of High grade engineering plastics, Colour Coded for easy identification, Calibration Conforming to Din 12650 & EN-ISO 8655. Standard ensuring high accuracy & precision. Does not require Calibration after After Autoclaving.



Cat No.	Volume μ l	Qty/pack	Price/Pc
F - 2.5	2.5	1	1500.00
F - 5	5	1	1500.00
F - 10	10	1	1500.00
F - 20	20	1	1500.00
F - 25	25	1	1500.00
F - 50	50	1	1500.00
F - 100	100	1	1500.00
F - 200	200	1	1500.00
F - 250	250	1	1500.00
F - 500	500	1	1500.00
F - 1000	1000	1	1500.00
F - 2000	2000	1	1500.00
F - 5000	5000	1	1500.00

LIQUID HANDLING PRODUCTS - BIOPETTE +

MICROPIPETTES, Fixed Volume Single channel, Fully Autoclavable, BIOPETTE Made of High grade engineering plastics, Colour Coded for easy identification, Calibration Conforming to Din 12650 & EN-ISO 8655. Standard ensuring high accuracy & precision. Does not require Calibration after After Autoclaving.



Cat No.	Volume μ l	Qty/pack	Price/Pc
V - 2.5	0.1 - 2.5	1	2750.00
V - 10	0.5 - 10	1	2750.00
V - 50	5 - 50	1	2750.00
V - 200	50 - 200	1	2750.00
V - 1000	200 - 1000	1	2750.00
V - 100	10 - 100	1	2750.00
V - 500	100 - 1000	1	2750.00
V - 20	2 - 20	1	2750.00
V - 300	20 - 200	1	2750.00
V - 5000	50 - 5000	1	2750.00
V - 10000	500 - 1000	1	2750.00

Pipette Controller - BIOPETTE +

Features :

- * Powerful pump with Double safety valve and Hydrophobic filter.
- * Soft grip push buttons.
- * Fill 25ml Pipette in just 4 seconds.
- * Safety valve and hydrophobic filter provide double protection against fluid penetration.
- * High and Low Speed, through varying Finger pressure.
- * Better Control of speed, through varying Finger pressure.
- * Low Battery Indicator.
- * Cadmium free environment friendly batteries.
- * Used Continuous for 8 Hours.
- * Fully Battery Charging time 2-3 hrs.
- * Pipette Volume Range from 0.1 - 200ml



Standard Accessories :

- * Battery Charger
- * Hydrophobic PTFE Membrane filter 0.45 μ l & 0.20 μ l each
- * Wall Mounting Stand with Screw.

Cat No.	Description	Range	Price/Pc
PC-01	Pipette Controller	0.1 - 200ml	8500.00



**P'fact Premium Pipette Tips Bulk 10µl Universal Graduated Tips
(1000 Tips / Pack)**

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-10B	-	-	-	500.00
PF-10LB	-	✓	-	820.00
PF-10FB	-	-	✓	1925.00
PF-10LFB	-	✓	✓	2300.00



**P'fact Premium Pipette Tips Bulk 10µl Long Universal Graduated Tips
(1000 Tips / Pack)**

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF -10XLB	-	-	-	500.00
PF -10XLLB	-	✓	-	820.00
PF -10XLFB	-	-	✓	1925.00
PF -10XLFB	-	✓	✓	2300.00



**P'fact Premium Pipette Tips Bulk 20µl Universal Graduated Tips (With Filter Only)
(1000 Tips / Pack)**

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-20FB	-	-	✓	1925.00
PF-20LFB	-	✓	✓	2300.00



**P'fact Premium Pipette Tips Bulk 50µl Universal Graduated Tips (With Filter Only)
(1000 Tips / Pack)**

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-50FB	-	-	✓	1925.00
PF-50LFB	-	✓	✓	2300.00



**P'fact Premium Pipette Tips Bulk 100µl Universal Graduated Tips (With Filter Only)
(1000 Tips / Pack)**

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-100FB	-	-	✓	1925.00
PF-100LFB	-	✓	✓	2300.00



**P'fact Premium Pipette Tips Bulk 200µl Universal Graduated Tips
(1000 Tips / Pack)**

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-200B	-	-	-	500.00
PF-200LB	-	✓	-	820.00
PF-200FB	-	-	✓	1925.00
PF-200LFB	-	✓	✓	2300.00



**P'fact Premium Pipette Tips Bulk 1000µl Universal Graduated Tips
(1000 Tips / Pack)**

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-1000B	-	-	-	550.00
PF-1000LB	-	✓	-	1010.00
PF-1000FB	-	-	✓	2290.00
PF-1000LFB	-	✓	✓	2700.00



P'fact Premium Pipette Tips Rack 10µl Universal Graduated Tips
(96 Tips x 10 Racks / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-10R	-	-	-	1690.00
PF-10SR	✓	-	-	1860.00
PF-10LR	-	✓	-	1960.00
PF-10FR	-	-	✓	2770.00
PF-10SLR	✓	✓	-	2130.00
PF-10SFR	✓	-	✓	2940.00
PF-10LFR	-	✓	✓	3000.00
PF-10SLFR	✓	✓	✓	3250.00



P'fact Premium Pipette Tips Rack 10µl Long Universal Graduated Tips
(96 Tips x 10 Racks / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-10LR	-	-	-	1690.00
PF-10LSR	✓	-	-	1860.00
PF-10LLR	-	✓	-	1960.00
PF-10LFR	-	-	✓	2770.00
PF-10LSLR	✓	✓	-	2130.00
PF-10LSFR	✓	-	✓	2940.00
PF-10LLFR	-	✓	✓	3000.00
PF-10LSLFR	✓	✓	✓	3250.00



P'fact Premium Pipette Tips Rack 20µl Long Universal Graduated Tips (with Filter only)
(96 Tips x 10 Racks / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-20FR	-	-	✓	2770.00
PF-20SFR	✓	-	✓	2940.00
PF-20LFR	-	✓	✓	3000.00
PF-20SLFR	✓	✓	✓	3250.00



P'fact Premium Pipette Tips Rack 50µl Long Universal Graduated Tips (with Filter only)
(96 Tips x 10 Racks / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-50FR	-	-	✓	2770.00
PF-50SFR	✓	-	✓	2940.00
PF-50LFR	-	✓	✓	3000.00
PF-50SLFR	✓	✓	✓	3250.00



P'fact Premium Pipette Tips Rack 100µl Long Universal Graduated Tips (with Filter only)
(96 Tips x 10 Racks / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-100FR	-	-	✓	2770.00
PF-100SFR	✓	-	✓	2940.00
PF-100LFR	-	✓	✓	3000.00
PF-100SLFR	✓	✓	✓	3250.00



P'fact Premium Pipette Tips Rack 200µl Universal Graduated Tips
(96 Tips x 10 Racks / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-200R	-	-	-	1690.00
PF-200SR	✓	-	-	1860.00
PF-200LR	-	✓	-	1960.00
PF-200FR	-	-	✓	2770.00
PF-200SLR	✓	✓	-	2130.00
PF-200SFR	✓	-	✓	2940.00
PF-200LFR	-	✓	✓	3000.00
PF-200SLFR	✓	✓	✓	3250.00



P'fact Premium Pipette Tips Rack 1000µl Long Universal Graduated Tips
(96 Tips x 6 Racks / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-1000R	-	-	-	1150.00
PF-1000SR	✓	-	-	1334.00
PF-1000LR	-	✓	-	1383.00
PF-1000FR	-	-	✓	1867.00
PF-1000SLR	✓	✓	-	1540.00
PF-1000SFR	✓	-	✓	2025.00
PF-1000LFR	-	✓	✓	2075.00
PF-1000SLFR	✓	✓	✓	2230.00



P'fact Premium Pipette Tips Reloading Stack 10µl Universal Graduated Tips
(96 Tips x 10 Plates / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-10K	-	-	-	1380.00
PF-10SK	✓	-	-	1500.00
PF-10LK	-	✓	-	1556.00
PF-10SLK	✓	✓	-	1790.00



P'fact Premium Pipette Tips Reloading Stack 200µl Universal Graduated Tips
(96 Tips x 10 Plates / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-200K	-	-	-	1380.00
PF-200SK	✓	-	-	1500.00
PF-200LK	-	✓	-	1556.00
PF-200SLK	✓	✓	-	1790.00



P'fact Premium Pipette Tips Reloading Stack 1000µl Universal Graduated Tips
(96 Tips x 5 Plates / Pack)

Cat.no	Sterile	Low Retention	Filter	Price/Pc
PF-1000K	-	-	-	805.00
PF-1000SK	✓	-	-	940.00
PF-1000LK	-	✓	-	980.00
PF-1000SLK	✓	✓	-	1120.00





P'fact Pipette Macro Volume Tips



Cat.no	Description	Packing	Qty/pack	Price/ Pc
PF-5000A	5000ul tip suitable to p'fact A & VA-FA Series pipette-2ml, 5ml, 500-5000ul & 1000-5000ul	Bulk	100 Tips	510.00
PF-10000A	10000ul tip suitable to p'fact A & VA-FA Series pipette-10ml. & 1000-10000ul pipette	Bulk	100 Tips	1500.00
PF-5000B	5000ul tip suitable to p'fact Prime 2 ml, 5ml & 500-5000ul pipette	Bulk	100 Tips	510.00
PF-10000B	10000ul tip suitable to p'fact Prime 10ml. & 1000-10000ul pipette	Bulk	100 Tips	1140.00



P'fact Tips Empty Racks

Cat.no	Description	Qty/pack	Price/ Pc
PF-10 -ER	Emptv Rack suitable for P'fact 10ul Tips	10 Racks	1125.00
PF-200 -ER	Emptv Rack suitable for P'fact 200ul Tips	10 Racks	1125.00
PF-1000 -ER	Emptv Rack suitable for P'fact 1000ul Tips	10 Racks	790.00

P'fact Centrifuge Tubes



Cat.no	Description	Packing	Price/ Pc
PCT-15-B	15 ml Centrifuge Tube with Cap	Bulk	3200.00
PCT-15-B-S	15 ml Centrifuge Tube with Cap.Sterile	Bulk	3700.00
PCT-15-R	15 ml Centrifuge Tube with Cap	Rack	4275.00
PCT-15-R-S	15 ml Centrifuge Tube with Cap.Sterile	Rack	4725.00
PCT-50-B	50 ml Centrifuge Tube with Cap	Bulk	4400.00
PCT-50-B-S	50 ml Centrifuge Tube with Cap.Sterile	Bulk	5300.00
PCT-50-R	50 ml Centrifuge Tube with Cap	Rack	6250.00
PCT-50-R-S	50 ml Centrifuge Tube with Cap.Sterile	Rack	6650.00
PCT-50S-B	50 ml Self Standing Centrifuge Tube with Cap	Bulk	4450.00
PCT-50S-B-S	50 ml Self Standing Centrifuge Tube with Cap.Sterile	Bulk	5400.00

P'fact Micro Centrifuge Tubes

Cat.no	Description	Packing	Qty/Pc	Price/ Pc
PCT-017B	1.7 ml Graduated Micro Centrifuge Tube	Bulk	1000 Tubes	650.00



INDEX

BRAKERS

Low Form	-	3
Tall Form	-	3

BOTTLES

B.O.D.	-	3
Laboratory Clear	-	4
Laboratory Amber	-	4
Reagent Narrow Mouth	-	3
Reagent Narrow Mouth Amber	-	4
Reagent Wide Mouth Amber	-	4
Reagent Wide Mouth Clear	-	4
Specific gravity. Glass Stopper	-	4
Specific gravity. Teflon Stopper	-	5
Weighing Tall form	-	5
Weighing Squat form	-	5

BURETTES

Automatic with Reservoir	-	6
Glass stop cock	-	5
Glass stop cock. schellbach	-	6
Rotoflow Teflon screwcock	-	5
Rotoflow Teflon Screwcock. schellbach	-	6
Teflon Stopcock	-	5
Teflon Stopcock. Schellbach	-	6

CONDENSERS

Air	-	6
Allihn	-	7
Double Surface	-	7
Graham	-	7
Liebig	-	6

CRUCIBLE, Filter

- 15

CYLINDERS

Measuring. Round base	-	7
Measuring. Hexagonal base	-	8
Measuring. Hexagonal base Schellbach	-	8
Nessler	-	7

DEAN & STARK APPARATUS

- 8

DESSICCATORS

Plain	-	8
Vacuum	-	8

DISHES, Petri

- 8

EXTRACTION APPARATUS

Soxhlet. Complete Set	-	9
Allihn Condenser. Spare	-	9
Extractor. Spare	-	9

FLASKS

Conical Narrow Mouth	-	9
Conical Wide Mouth	-	9
Conical Wide I/C Stopper	-	9
Distilling	-	12
Filter	-	12
Flat Bottom Plain	-	10
Flat Bottom. Centre neck with I/C Joint	-	10-11
Iodine	-	10
Kjeldahl. Plain	-	13
Kjeldahl. with I/c Joint	-	13
Round Bottom Plain	-	10
Round Bottom Two Necks. Annealed	-	11
Round Bottom Two Necks. Parallel	-	12
Volumetric. Clear	-	13
Volumetric. Amber	-	13

FUNNELS

Buchner	-	15
Long Stem	-	13
Separating. Pear Rotoflow Teflon Screwcock	-	14
Separating. Pear Glass Stopcock	-	14
Separating. Pear Teflon Stopcock	-	14
Separating. Globe. Rotoflow Teflon Screwcock	-	14
Separating. Globe. Glass Stopcock	-	14
Separating. Globe. Teflon Stopcock	-	14
Short Stem	-	13

LIQUID HANDLING PRODUCTS

Micropipettes. Fixed Volume BIOPETTE	-	18
Micropipettes. Variable Volume BIOPETTE	-	18
Pipette Controller BIOPETTE	-	18

MEMBRANE FILTER ASSEMBLY


- 17

PIPETTES

Mohr's	-	15
Serological	-	15
Volumetric	-	16

Terms of Sale

- ❖ Prices Mentioned in the Price List are the Maximum recommended Prices per Unit/ Piece in Indian Rupees.
 - ❖ Prices are Ex- Godown, Chennai, Packing, Forwarding Expenses extra at Cost. VAT or CST Extra as per Govt. rules.
 - ❖ Prices are Subject to change without Notice. Prices Prevailing at the time of actual delivery will be Valid.
 - ❖ Invoices are to be paid as per agreed Credit period without any deductions.
In Case of outstation Cheque / Drafts, The Bank Charges will be Customers account.
 - ❖ The Prices are for the Standard Pack for each item, it would be appreciated that Customer places orders for the items based on the Standard pack Size indicated and in multiples thereof.
 - ❖ No goods should be returned to us without our Prior approval.
 - ❖ The Sales of the laboratory glassware products will be Subject to Chennai Jurisdiction for all the Parties and for all legal matters, the Indian Laws are applicable.
-
-





CHENNAI

7A, Natarajan St., (Saidai Sambadhan St.)
Saidapet, Chennai-15.
Phone : 24714329, Tele Fax : 044-43302975
Cell : 98400 54224, 99400 94224
Email : kavinscientific@yahoo.co.uk

MADURAI

233, South Marat Street,
South Gate, Madurai-625001
Cell : 98404 54224
Email : kavinscientific.mdu@gmail.com



AUTHORISED STOCKIST :

